





# Acry Lux V

# Acry Lux



# Anterior Upper



Combined Articulations	Anterior			Posterior	
	Upper	Lower	Lower		
E1	K5	K3		M2	
E2	K6	K5		M4	N3
E3	K6			M4	N3
E4	K9			M5	
E5	K10			M5	
F1	K1			M1	N1
F2	K2			M2	
F3	K5	K2		M3	N2
F4	K10	K9		M5	
G1	K2	K1		M1	N1
G2	K6	K7		M2	
G3	K6	K4		M4	N3
G4	K8	K1		M5	
H1	K7	K1		M1	N1
H2	K5	K2		M2	
H3	K3	K1		M3	N2
H4	K4	K3	K5	M4	N3

Combined Articulations	Anterior			Posterior	
	Upper	Lower	Lower		
I1	K7	K1		M1	N1
I2	K5			M3	N2
I3	K11			M4	N3
I4	K8	K12		M5	
L1	K8			M4	N3
L2	K11			M4	N3
L3	K11			M4	N3
L4	K12			M5	
32	K3	K1		M3	
52	K6			M3	M4
62	K6	K4		M3	
66	K6	K5		M3	M2
73	K10	K9		M5	M4

Due to lighting conditions during photographic process the product illustrated may vary in shade and mould comparison. For shade and mould verification please refer to the live guides.

# Anterior Lower

I1	1	42.3
	2	9.3
	3	11.5
	4	7.9



I2	1	44.3
	2	9.7
	3	11.5
	4	8.0



I3	1	45.2
	2	9.0
	3	11.4
	4	8.3



I4	1	47.0
	2	9.7
	3	11.5
	4	8.7



L1	1	47.0
	2	9.3
	3	11.2
	4	8.2



L2	1	47.1
	2	9.3
	3	11.7
	4	8.7



L3	1	48.1
	2	9.7
	3	11.4
	4	9.0



L4	1	50.5
	2	9.9
	3	12.0
	4	9.2



32	1	40.4
	2	8.7
	3	10.5
	4	7.7



52	1	43.7
	2	10.4
	3	12.5
	4	8.2



62	1	47.1
	2	9.7
	3	11.5
	4	8.8




66	1	47.7
	2	10.6
	3	12.0
	4	8.4



73	1	48.0
	2	11.9
	3	14.0
	4	8.7



K1	1	32.3
	2	8.1
	3	9.7
	4	4.4



K2	1	32.7
	2	10.0
	3	11.4
	4	4.5



K3	1	33.8
	2	8.5
	3	9.7
	4	4.6



K4	1	34.6
	2	9.8
	3	11.3
	4	5.0



K5	1	34.4
	2	9.0
	3	10.9
	4	5.1



K6	1	35.0
	2	9.8
	3	11.0
	4	5.3



K7	1	34.8
	2	8.4
	3	9.6
	4	5.0



K8	1	40.0
	2	8.2
	3	10.1
	4	6.1



K9	1	38.5
	2	9.5
	3	11.2
	4	6.5




K10	1	40.0
	2	12.2
	3	14.0
	4	5.7




K11	1	39.4
	2	8.8
	3	10.8
	4	5.8



K12	1	41.2
	2	9.4
	3	11.4
	4	6.2

# Posterior Upper - Lower



# Acry Lux

# Acry Lux V

# Acry Lux

- Three layer acrylic tooth
- 19 Ivoclar PE comparative shades
- Available in:
  - 30 upper anterior moulds
  - 12 lower anterior moulds
  - 8 posterior moulds
 (Including variation in cusp steepness)


# Acry Lux V

- Three layer acrylic tooth
- 15 Vita Lumin Vacuum® comparative shades
- Available in:
  - 30 upper anterior moulds
  - 12 lower anterior moulds
  - 8 posterior moulds
 (Including variation in cusp steepness)

Acry Lux		Comparative	
Shade	Colour	Shade	Brand
1A	White		
2A	White		
3A	White		
4A	Grey		
2B	Yellow		
4B	Grey		
5B	Dark Yellow		
6B	Grey		
1C	Yellow		
2C	Yellow		
3C	Brown		
4C	Brown		
6C	Grey		
1D	Brown		
4D	Brown		
6D	Brown		
1E	Yellow		
2E	Dark Yellow		
3E	Dark Yellow		


Acry Lux V		Comparative	
Shade	Colour	Shade	Brand
A1	White		
A2	Brown		
A3	Brown		
A3.5	Brown		
A4	Brown		
B2	Yellow		
B3	Yellow		
B4	Yellow		
C1	Grey		
C2	Grey		
C3	Grey		
C4	Grey		
D2	Reddish Grey		
D3	Reddish Grey		
D4	Reddish Grey		

## Mould Specifications



**Anteriors Upper - Lower**

- > All specifications in millimetres (mm)
- 1 Width of anterior set
- 2 Length of central incisor excluding root surface
- 3 Length of central incisor including root surface
- 4 Width of central incisor



**Posterior Upper - Lower**

- > All specifications in millimetres (mm)
- 1 Width of upper posterior set
- 2 Width of lower posterior set